

Experiment - I

Aim - To understand the phenomenon Photoelectric effect and plot a graph between photocurrent & potential to determine stopping potential.

Theory -

The photoelectric effect is an instantaneous phenomenon. There is no time delay b/w the incidence of light and emission of photo electrons.

The no. of photoelectrons emitted is proportional to the q incident light. Also, energy of ~~incident light~~ photoelectron is independent of incident light.

The energy of emitted photoelectrons is directly proportional to the frequency of light.

It has been observed that there must be a minimum energy needed for electrons to escape from a particular metal surface and is called work function ' W '.

$$W = h\nu_0$$


acc. to Einstein Photoelectric effect equation

$$h\nu = KE + W \quad -$$

$$KE = h\nu - \nu_0$$

If we increase the reverse potential, the photocurrent gradually decreases and becomes zero, this is the minⁿ applied reverse potential called stopping potential V_0 .

$$KE_{max} = eV_0$$



Set up for
Photoelectric Effect

Procedure -

- 1 Select the material for studying photoelectric current
- 2 Select the area of the material, wavelength, intensity of incident light.
- 3 Switch on the light source
- 4 Measure the reverse current for various reverse voltage
- 5 Plot the current-voltage graph and determine the threshold voltage.
- 6 Repeat the experiment by varying the intensity for a particular wavelength of incident light.
- 7 Repeat the experiment by varying the wavelength for a particular intensity of the incident light.


Result -

Metal used - Copper

~~Stop~~

Threshold ~~voltage~~ frequency = $1.1 \times 10^{15} \text{ Hz}$

Stopping potential = -7.8 V



22/6/22

Observations -

constant intensity ($I = 5 \text{ W/m}^2$) different wavelength

Cu

$\lambda = 100\text{nm}$		$\lambda = 150\text{nm}$		$\lambda = 200\text{nm}$		$I_p (\mu\text{A})$
V	$I_p (\mu\text{A})$	V	$I_p (\mu\text{A})$	V	$I_p (\mu\text{A})$	
-7.8	0	-3.6	0	-1.6	0	
-7	1.79	-3	1.44	-1	1.21	
-6	4.29	-2.5	2.69	-0.5	2.52	
-5	6.79	-2	3.94	0	3.77	
-4	9.29	-1.5	5.19			
-3	11.79	-1	6.49			
-2	14.29	0	8.94			
-1	16.79					
0	19.29					

for constant wavelength $\lambda = 100\text{nm}$, variable intensity



Intensity = 5 W/m^2	Intensity = 10 W/m^2	Intensity = 15 W/m^2
0	19.29	0
-1	16.79	-1
-2	14.29	-2
-3	11.79	-3
-4	9.29	-4
-5	6.79	-5
-6	4.29	-6
-7	1.79	-7
-7.8	0	-7.8

$$\text{Calculations} - E = h\nu + kT \rightarrow h\nu = \omega_0 + \epsilon kT$$

$$\omega_0 = \frac{6.6 \times 10^{-34} \times 5 \times 10^3}{10^7} - (1.6 \times 10^{-19})(-7.8)$$

$$\omega_0 = 7.32 \times 10^{-19} \text{ J} = 4.57 \text{ eV}$$

$$\nu_0 = \frac{7.32 \times 10^{-19}}{6.6 \times 10^{-34}} = 1.1 \times 10^{15} \text{ Hz}$$



Aim- To determine the specific heat of Copper, Lead and glass.

Materials Required-

Dewar flask, steam generator, heating chamber, beaker, metal samples, silicon tube,

Theory-

Heat capacity C of an object is defined as $C = \frac{\Delta Q}{\Delta T}$, where ΔQ is the amount of heat required to change the temperature of object by ΔT . Specific heat is the heat capacity per unit mass.

$$C = \frac{c}{m} = \frac{\Delta Q}{m \Delta T}$$

$$\Delta Q = Cm \Delta T$$

According to conservation of energy -

$$\Delta Q_{\text{gained}} = \Delta Q_{\text{lost}}$$

$$m_1 c_1 (T_f - T_i) = m_2 c_2 (T_f - T_i)$$

for this experiment, the system consist of mixing a given mass m_1 of a hot metal with specific heat c_1 at temperature T_i and a known mass m_2 of water with specific heat c_2 at a lower temperature T_i , contained in calorimeter of mass m_3 , c_3 - specific heat of calorimeter, initial temperature T_i . Let final temp be T_f .

$$\Delta Q_{\text{lost}} (\text{metal}) = \Delta Q_{\text{gained}} (\text{water}) + \Delta Q_{\text{gained}} (\text{calorimeter})$$

$$m_1 c_1 (T_i - T_f) = m_2 c_2 (T_f - T_i) + m_3 c_3 (T_f - T_i)$$

Teacher's Signature :

$$m_p c_1 (t_s - t_b) = (M_w + W) (t_b - t_2)$$

W = water equivalent - M_w is the mass of water which absorb the same amount of heat as absorbed by flask to raise its temp by 1°C .

We will consider a metal specimen heated to a high temperature is dropped into water contained in a calorimeter at a lower temperature. Specific heat can be obtained by equating the heat lost by the metal raised by both calorimeter and water.

Procedure -

- 1 The experiment setup is placed
- 2 Open the top of steam generator by turning it anticlockwise and add sufficient water so that heater element is well dipped into water.
- 3 Connect the heating chamber to the steam generator using silicon tube.
- 4 Attach silicon tubing to the bottom base connection of the heating chamber and hang the other end in the beaker. See that the silicon tubing is securely sealed at all connections.
- 5 Fill the sample chamber of the heating chamber with a weighed quantity of lead shots and seal it with stopper.
- 6 Turn on the steam generator to generate steam and heat the shots for about 20-30 minutes in the heating chamber.

Determination of specific heat

- 1 Empty the dewar flask and place it in water so that it regain its normal temperature.
- 2 Open the cover of the dewar vessel and shift below the heat chamber and drop the shots at 100°C into the dewar flask.
- 3 Stir the mixture gently by rotating the flask clockwise and anti-clockwise with hand.

Observation -

Determination of water equivalent of dewar flask -

1. mass of flask (M_1) =
2. M_2 (mass of cold water (M_2) + mass of flask (M_1)) =
3. mass of cold water (M_2) = $M_2 - M_1$ =
4. Temp of cold water (T_c) =
5. Temp of hot water (T_h) =
6. M_3 (cold water + flask + hot water) =
7. Mass of hot water = M_h = $M_3 - M_1 - M_2$ =
8. Temp of mixture = T_m =

$$\text{heat lost} = \text{heat gained}$$

$$M_h (T_h - T_m) = M_c (T_m - T_c) + w (T_m - T_c)$$

$$w =$$

Determination of the specific heat -

- (i) Sample result for lead -
- initial temperature of water (T_2) =
- final temperature of water (T_b) =
- mass of lead shot, m_{pb} =
- mass of water in flask in water with lead shot) =
- water equivalent =
- t_s = temp of steam =

$$m_{pb} C_1 (t_s - t_1) = (M_w + w) (T_b - T_2)$$

Literature value = $0.129 \text{ kJ kg}^{-1} \text{ }^{\circ}\text{C}^{-1}$

$$\text{expt} =$$

- 4 Read the temperature of the mixture when the temperature of the water stop rising and use eqn - 6 to find specific heat.
5 Rearrange the experiment with copper and glass slab.

Result

specific heat of lead =

specific heat of copper =

specific heat of brass =

(ii) Sample result for Copper -
initial temp of water (t_2) =
final temp of water (C_4P) =
 m_w = mass of Cu Sheet, M_w =
 M_w = mass of water + Cu =
Water equivalent W =
 k_s = temp of steam

$$m_w C_4(t_5 - t_4) = (M_w + W)(C_4 - t_2)$$

C_2 =
literature value =

error =

(iii) Sample result for glass -

initial temp of water (t_2) =
final temp of water (C_4P) =
mass of glass block, M_w =

M_w =

w =

k_s =

$$m_w C_3 (t_5 - t_4) = (M_w - w)(C_4 - t_2)$$

C_3 =
literature value =
error =

Aim -

To determine the wavelength of a laser using the Michelson interferometer

Apparatus-
laser light source, Michelson interferometer kit, optical bench, meter scale

Theory-

In Michelson interferometer, colored beams are obtained by splitting a beam ~~because~~ of light using beam splitter. The rays falling on mirror M_1 and M_2 reflect and interfere on screen. Two virtual images are formed. S_1 due to reflection at M_1 , and S_2 due to reflection at M_2 separated by a distance

2d


Using Michelson interferometer, wavelength of light λ is determined by

$$\lambda = \frac{2d}{N}$$

Procedure-

- 1 The laser beam must strike at the center of movable mirror and should be reflected directly back into the laser aperture
- 2 Adjust the position of beam splitter so that the beam is reflected the fixed mirror.
- 3 Adjust the angle of the beam splitter to make the two sets of spots as close together as possible.
- 4 With the screen on the back of the adjustable mirror adjust the mirror's tilt until the two sets of spots coincide
- 5 Expand the laser beam slowly rotating the collimating lens in front of laser

Teacher's Signature: _____



- 5 Align the laser with the micrometer and make sure that the fingers are moving when the micrometer screw is turned.
- 7 Make a point on the screen and note the micrometer reading
- 8 Calculate the distance d , repeat the procedure to take average distance.
- 9 With a known wave length, $\lambda = \frac{c}{n\Delta}$ is calculated.
- 10 Unknown wavelength is calculated.

Result -
Second
The wavelength of the given ~~laser~~^{light source} is 550 nm.

✓
27/6/22

Observation-

Least count = 0.001 mm

L.C.CSP = $\frac{0.01}{\sqrt{N}}$

Distance measured for N trials d = 10 ✓

No.	CSL (mm)	VSL (mm)	Distance d (mm)	Difference (mm)
1	90×0.0001	0.5×0.1	0.0140	0.0031
2	59×0.0001	0.5×0.1	0.0109	0.0029
3	30×0.0001	0.5×0.1	0.0080	0.0028
4	12×0.0001	0.4×0.1	0.0052	0.0022
5	0×0.0001	0.5×0.1	0.0030	

Calculation-

$$\text{mean} = \Delta d = \frac{0.0031 + 0.0029 + 0.0028 + 0.0022}{4}$$

$$= 0.00275$$

$$\lambda = \frac{2d}{N} = \frac{2(0.00275)}{10} = 0.00055 \text{ mm}$$

$$= 550 \text{ nm}$$

$$\text{Error} = \frac{620 - 550}{620} \times 100 = 11.1\% \cancel{\text{error}}$$

Aim-

To determine the wavelength of laser light using diffraction grating.

Theory- Diffraction of light occurs as it passes through a slit, the angle to the minima in the diffraction pattern is given by -

$$a \sin \theta = m\lambda$$

where, a = slit width

θ = angle from the centre of pattern to the m^{th} minima.

λ = wavelength of the light and

m = order.

Since, the angles are usually small, it can be assumed that

$$\tan \theta = \frac{y}{D}$$

Where, y = distance on the screen from the centre of the pattern to the m^{th} minimum

D = distance from the slit to the screen

$$\therefore \text{Slit width } a = \frac{m \lambda D}{y} \quad (m = 1, 2, 3, \dots)$$


Procedure 1 Set up the diode laser at the end of the optical bench and place the slit holder about 5 cm in front of laser.

2 Put the photo-detector at the other end of the ~~optical~~ optical bench adjust the vertical and horizontal position and laser beam so that beam falls at the centre of the pin hole.

3 Connect the ~~multimeter~~ multimeter to the detectors.

4 Mount the single slit in slit holder and adjust the position of laser beam from left to right and up and down until the beam is centred on the phot.

Teacher's Signature : _____



Figur. Diffraction grating setup.

- 5 Determine the distance b/w slit and detector
6 Record the intensity at closed intervals by traversing the detector through the cross-section.

- 7 Diode laser is mounted on the Saddle
8 A plane transmitting grating is mounted on an upright next to laser.
9 The position of α of the spot of 1st order on either side of central maxima is marked
10 The distance D b/w the grating and screen is measured
11 Calculate $\sin\theta = \frac{\alpha}{D}$

12 Grating element d for grating lines per mm
 $= \frac{1}{\lambda \times 100}$

Result -

The wavelength of the given laser beam = 642.75 nm

~~Ques 12~~

Observations

$$D = 31 \text{ cm}$$

Sno.	Diameter (mm)	Left (mm)	Right (mm)	R (mm)	D (mm)
		$N = 300 \text{ lines/mm}$			
1	1	10.55	10.5	0.025	11
2	2	7.74	7.81	0.035	9
		$N = 80 \text{ lines/mm}$			
3	1	10.4	10.2	0.2	23.7
4	2	24.615	25.20	0.375	23.7

卷之二

$$\lambda_1 = \frac{300(0.025)}{(1)(1)} = 681 \text{ nm}$$

$$g_2 = \frac{300 \times 0.032}{(2)(\omega)} = 583 \text{ mm}$$

$$\lambda_3 = \frac{(\text{波長} 0.2)}{(1)(23.7)} = 615 \text{ nm}$$

$$\lambda_4 = \frac{(80)(0.379)}{(2)(23.7)} = 6.32 \text{ nm}$$

$$\underline{\text{Mean}} \quad \lambda = \frac{\lambda_1 + \lambda_2 + \lambda_3 + \lambda_4}{4} = 642.75 \text{ mm}$$

Calculations

~~22/02/2022~~

Aim -

To determine the resistivity of semiconductors by four probe method.

Apparatus -

four probe arrangement, Sample, oven 0-200°C, constant current generator, oven power supply, digital panel meter.

Theory -

$$\text{Resistance} \quad R = \rho \frac{L}{A} \quad \text{where } \rho \text{ is resistivity}$$

~~Resistivity of material is uniform in the area of measured if there is no minority carrier injection into the semiconductor by the current carrying electrode, they recombine so that effect on conductivity is low.~~

A semiconductor has electrical conductivity intermediate in magnitude. In case of intrinsic semiconductors, the Fermi level lies in between conduction and valence band. If resistance is infinite. As temperature increases, the occupancy of ~~resistor~~ conduction band goes up and results in decrease in electrical resistivity of semiconductors.

$$\rho = \frac{\rho_0}{6(\omega/s)} \quad \rho_0 = \frac{V}{I} \quad 2\pi s \quad S = \text{spacings between probes}$$

Procedure -

- 1 Select the Semiconductor material from the combobox.
- 2 Select the Source current from the slider. Restrict the slider based on current range.
- 3 Select the range of oven from the combobox.
- 4 Set the temperature from the slider.

- 5 Click on the Run button to start heating the oven in a particular interval, from the default 28°C to the temp we set and then stop heating when temperature is reached.
- 6 Measure to display the temp in the oven.
- 7 Select range of voltmeter and measure voltage.
- 8 Calculate the resistivity of semiconductor in eV for given temperature.
- 9 Temperature vs Resistivity graph is plotted

Result -

The resistivity of the given semiconductor = 6.204 Ohm cm

~~12/16/22~~

Observations -

Temperature $T, (K)$	Voltage $V, (mV)$	Current $I, (mA)$	Resistivity ρ (ohm cm)
25	87.24	3	6.204
30	84.65	3	6.020
35	82.22	3	5.845
40	79.04	3	5.683
45	77.78	3	5.530
50	75.75	3	5.385
55	73.83	3	5.249
60	72.01	3	5.119
65	70.30	3	4.998
70	68.67	3	4.882


Calculations

Distance between probes ; $s = 0.2 \text{ cm}$
 Thickness of sample ; $w = 0.05 \text{ cm}$
 from standard table $b(w_{15}) = 5.89$

$$\rho = \frac{V}{I} = 36.54 \Omega \text{ cm}$$

$$P = \frac{\rho}{l} = 1.85 \Omega \text{ cm}$$

Graph -



Aim -

To find the numerical aperture of a given optic fibre and hence to find its acceptance angle

Materials required

Optical fibre cable, cutter, output unit, screen,

Theory-

Numerical aperture of any optical system is a measure of how much light can be collected by optical system. It is product of refractive index of incident medium and sine of maximum ray angle.

$$NA = n \sin \theta_{\max}; \text{ but air is 1}$$

$$\therefore NA = \sin \theta_{\max}$$

for a step index fibre, as in present case, numerical aperture is given by

$$NA = \sqrt{n_{core}^2 - n_{cladding}^2}$$

for very small difference in refractive indices, equation reduces to

$$NA = n_{core} \sqrt{2\Delta}$$

where Δ is fractional difference in refractive indices.

The fibre may refer to specification of Park fibre and record many values NA, n_{core}, n_{cladding}.

Procedure-

- 1 Connect one end of cable 1 to ~~F0~~ LED in fibre optic LED driver section of FO & and other end to N.A. fig.
 - 2 Switch on driver. Light should appear at end of fibre or N.A. fig. Turn set Pot knob clockwise to set to maximum po. The light intensity should be increased.
 - 3 Hold white screen with ~~concave~~ concave concentric circles vertically at a suitable distance to make red spot from emitting fibre coincide with 10 mm circle. Note that converging circle must coincide with light. Record 'L' distance of green from fibre end and note diameter of spot.
 - 4 Compute N.A. from formula—
- $$\text{N.A.} = \sin \theta_{max} = \frac{w}{\sqrt{4L^2 + w^2}}$$
- 5 In case fibre is under filled, intensity within spot may not be evenly distribution. To ensure even distribution of light in fibre, first remove twists on fibre and then wind 5 turns of fibre on to mandrel.

Result-

Numerical aperture of green wire is = 0.69 radian

Observation -

S. no.	L (cm) (approx)	W (cm) diameter	sin θ (N.A)
1	1	1.7	0.640
2	2	2.9	0.587
3	3	6.4	0.530
4	5	7.8	0.870
5	8	10.4	0.460
6	10		

$$\text{mean } \bar{\theta} = N \cdot N = 0.61$$

$$N = A = 0.61$$

$$\theta_{\text{max}} = 0.69 \text{ radian}$$

Aim- To determine wavelength of sodium light by
Newton rings

Apparatus required- Plane-concave lens of large radius of curvature, optical arrangement for Newton's rings, plane glass plate, sodium lamp, and travelling microscope.

Theory-

Due to air film formed by glass plate & the plane concave lens of large radius of curvature, interference fringes are formed which are observed directly through travelling microscope. Rings are concentric circles path difference : $\Delta n = 2 \pi R \Delta \phi$

~~for wavelength~~


$$\text{using formula: } \lambda = \frac{D^2 m + n - D_n^2}{4 n R}$$

- $D_{m+n} \Rightarrow$ diameter of $(m+n)^{\text{th}}$ ~~ring~~ ring
- $D_n \Rightarrow$ diameter of n^{th} ring
- ~~m~~ \Rightarrow an integer of no. of rings
- $R \Rightarrow$ ~~Radius~~ curvature of plane - convex lens
- $n \Rightarrow$ order of ring

Result

Procedure-

- 1 After the experimental arrangement, the glass plate is inclined at an angle 45° to the horizontal. The glass plate reflects light from the source vertically downwards and falls normally on the convex lens.



apparatus for Newton rings

2. Michelson rings are seen using a long focused microscope, focussed on the air film. The ~~co~~ vernier microscope is made tangential to the 1st ring on the left side of the centre. The reading of main scale and vernier scale are noted.
3. Then it is seen through liquid side from the centre and observation are taken. Similarly readings for 2nd, 3rd, 4th etc. ring are taken.
4. Diameter of ring is found out by subtracting readings on the left and right sides. The squares of diameter and hence $D^2 - D_{\text{new}}$ are found out. Wavelength for individual rings were calculated and their mean was found.

Result -

~~Wavelength of Sodium Lamp is = 542 nm~~

~~Actual calculated value of wavelength = 589 nm~~

Percentage error = 7.9%

Ram

Observations -

R = 60 mm

Order of avg (P)	Micrometer reading		Diameter D _P (mm) left-right	D _P ² (mm ²)	D _{mp} ² - D _P ² (mm ²) m=2	$\frac{2}{D_{mp}^2 + D_p^2}$ LPR
	left	right				
1	10.365	10.60	-0.235	0.056225	0.01775	369
2	10.60	10.24	0.36	0.1296	6.994	207
3	10.315	10.06	0.755	0.57	1.3621	283
4	10.99	9.93	1.06	1.1236	14.846	309
5	11.14	9.45	1.34	1.921	6.1335	1277
6	11.23	9.65	1.615	2.608225	6.3617	1325
7	11.33	8.69	2.84	8.0656	1.44433	300
8	11.15	8.42	2.945	8.97	1.9801	413
9	11.49	8.32	3.17	10.0489	2.1312	444
10	11.55	8.24	3.31	10.4661	2.364	493
11	11.65	8.16	3.49	12.1801		
12	11.76	8.11	3.65	13.3225		

Calculation

$$\text{Mean} = \frac{S_{420}}{10} = 542 \text{ mm}$$

$$\text{Error} \Rightarrow \frac{S_{89} - S_{42}}{S_{89}} \times 100 = 7.9\%$$



Experiment -8

Aim- To determine the value of g , the acceleration due to gravity at a particular location using Kater pendulum.

Apparatus required -
Kater's pendulum, stopwatch, meter scale, weights, wooden and metal souffle cones.

Theory -

We consider the force of gravity to be acting at G. If h_i is the distance to a from the suspension point O_i at the knife edge K_i

$$T_i \ddot{\theta} = -mg h_i \sin \theta$$

where T_i is moment of inertia.

Comparing to the motion of simple pendulum

$$M T_i^2 \ddot{\theta} = -mg h_i \sin \theta$$

$$M g h_i / T_i = -g / \omega$$

$$\text{or } T_i = M k_i^2 \quad \text{--- (1)}$$

$$k_i^2 = h_i l_i \cdot \text{--- (2)}$$

If I_g is the moment of inertia of pendulum about G, radius of gyration is $I_g / M k_i^2$


$$l_i = \frac{l_1^2 + l_2^2}{h_i} \quad \text{using (2)}$$

$$\text{true period } T_i = 2\pi \sqrt{\frac{l_i}{g}} = 2\pi \sqrt{\frac{I_g^2 + k_i^2}{g h_i}} \quad \text{--- (3)}$$

$$\text{Squaring (3) and calculate } g = \frac{8\pi^2}{h_1 + h_2} \left[\frac{l_1^2}{h_1^2} + \frac{l_2^2}{h_2^2} \right]^{-1}$$

Teacher's signature : _____

Experimental Setup for Kotsu Pendulum



Procedure -

1. Choose a desired answer from pendulum.
2. Shift the weight w , to one end of Kater's pendulum and fix it. Fix the knife edge k_1 just below it. Keep the knife edge k_2 at the other end. Keep the small weight w near to centre.

~~W~~ suspend the pendulum about the knife edge and take the time for about 10 oscillations. Note the time t_1 using a stopwatch and calculate its time period using $T_1 = t_1 / 10$.

3. Now suspend about knife edge k_2 by inverting the pendulum and note the time t_2 for oscillation. Calculate T_2 .
4. If $T_2 \neq T_1$, adjust the position of knife edge k_2 so that $T_2 = T_1$. Balance the pendulum on a sharp edge and mark the position of its centre of gravity. Measure the distance of the knife edge k_2 and that of k_2 as l_2 from k_1 .

Result

The accn due to gravity at a given place is found

$$\text{Accn} = 10.38 \text{ m s}^{-2}$$

$$\text{Percentage error} = 5.91\%$$

Observation -

S.no	Distance (cm)	T ₁ Period	T ₁	T ₂ Period	T ₂	g
1	79.5	17.28	1.728	17.48	1.748	10.38

Calculations -

T₁ = T₂

$$g = \frac{8\pi^2 (L_1 + L_2)}{C(T_1^2 + T_2^2)}$$

$$= \frac{8 \times (3.14)^2 \times 79.5}{2.988 + 3.055}$$

$$= \frac{6270.7}{6.040}$$

$$= 10.38 \text{ cm/s}^2 = 10.38 \text{ m/s}^2$$

$$\% \text{ error} = \frac{10.38 - 9.8}{9.8} \times 100 = 5.91 \%$$

Aim -
To measure the height of window using Sextant

Materials Required
Sextant, measuring tape, Rigid clamp stand

Theory -

height of accessible object -

Let h = height

θ = angle made by PQ at O

$$\therefore h = d \tan \theta$$

d = OP measured using measuring tape

Determination of height of inaccessible object -

measure angle θ_1 and θ_2 at R and S using Sextant

$$h = PQ =$$

$$\therefore \cot \theta_1 = \frac{PR}{h} \quad \cot \theta_2 = \frac{SP}{h}$$

$$\cot \theta_2 - \cot \theta_1 = \frac{SP}{h} = \frac{d}{h}$$




$$h = d / (\cot \theta_2 - \cot \theta_1)$$

Determination of height of a tall object -

$$\text{height of object} = QP - SP$$

$$= d(\tan \theta_2 - \tan \theta_1)$$

$$\therefore h = d(\tan \theta_2 - \tan \theta_1) \quad \text{①}$$



Procedure-

- 1 Find the vernier constant of scale
- 2 In the line of your eye, mark a horizontal line on the tower. Place the sextant at a large distance from source.
- 3 Observe the horizontal line and observe its one image directly and other by two reflections at M_1 & M_2 . Now rotate the movable arm, such that the horizontal line is seen continuously. If there is no zero error, zero reading should coincide with zero of vernier scale. If not, record the zero error with protractor sign.
- 4 Now turn the movable arm slowly till the image of top window with the image of line. Record this reading also. The difference of two readings is the angular elevation θ_1 .
- 5 From this point of observation move the sextant to about 3-4.5 m away.
- 6 Repeat the step 4 & 5 find angle θ_2 .
- 7 Using eqn 1 to find the height of the window. Also measure the distance of horizontal line from the ground. Add this in the calculation of actual height.

Result-

height of window calculated using measuring tape = 135 cm
 height of window calculated using sextant = 139.24 cm
 percentage error = 3%

Teacher's Signature : _____

Observation-

$$D = 387.35 \text{ cm}$$

$$\text{Lc of main scale} = 1^{\circ}$$

$$\text{Lc of circular scale} = \left(\frac{1}{60}\right)^{\circ} = 1'$$

Vernier Constant Calculation:

$$4 \text{ CSD} = 5 \text{ VSD}$$

$$1 \text{ VSD} = \frac{4}{5} \text{ CSD} = \frac{4}{5} \left(\frac{1}{60}\right)^{\circ} =$$

$$VC = 1 - VSD = \left(1 - \frac{4}{5}\right) \left(\frac{1}{60}\right)^{\circ} = \left(\frac{1}{5}\right) \left(\frac{1}{60}\right)^{\circ} = \left(\frac{1}{300}\right)^{\circ} = 12'$$

Log of Vernier scale = $(1/300)^{\circ}$

S.No	MSR	CSR	VSR	\oplus (deg)
1	14	$15 \times \frac{1}{60} = 0.25$	$4 \times \frac{1}{300} = 0.0133$	14.2633
2	31	$27 \times \frac{1}{60} = 0.45$	$5 \times \frac{1}{300} = 0.0166$	31.4666

Calculations

$$\begin{aligned} h &= D (\tan \theta_2 - \tan \theta_1) \\ &= 387.35 (\tan (31.46) - \tan (14.26)) \\ &\approx 139.24 \text{ cm} \end{aligned}$$

$$\text{Error} \rightarrow \frac{139.24 - 135}{135} \times 100 = 3\%$$